

Invacare® Respiratory Products



Yes, you can.™



Invacare® Piston Regulators

Invacare offers the right selection of regulator models, sizes and styles. Plus, the 6-year warranty ensures years of accurate, trouble-free service.

FEATURES

- “Select flow” style
- Lightweight and durable
- Bumper-protected contents gauge
- Flow lock knob
- Competitively priced

Model Selection Guide

| Model no. | Units/Case | Warranty | Description | Maximum flow (lpm) |
|-----------|------------|----------|--|--------------------|
| IRC2121 | 1 | 6 years | CGA870, 0-6 lpm | 6 |
| IRC2020 | 1 | 6 years | CGA870, 0-8 lpm | 8 |
| IRC2025 | 1 | 6 years | CGA870, 0-15 lpm | 15 |
| IRC2002 | 1 | 6 years | CGA870, 0-4 lpm (in pediatric increments) | 4 |
| IRC2030 | 1 | 6 years | CGA540, 0-15 lpm | 15 |



Invacare® Compressed Oxygen Systems

Invacare compressed oxygen systems meet or exceed the requirements to support the home oxygen therapy patient and provider. They are applicable for virtually all home oxygen prescriptions. And because they contain fewer parts, Invacare systems minimize inventory space.

FEATURES

- Lightweight aluminum cylinder with contents gauge on regulator
- Two-wheel cart or black nylon shoulder bag
- Durable 0-8 lpm Invacare piston regulator to match cylinder
- Nasal cannula
- Plastic wrench
- Operating instructions
- Lightweight shoulder pack weighs just 8 lb.

Model Selection Guide

| Model no. | Units/Case | Description |
|-----------|------------|--|
| HP4121 | 1 | D-size shoulder pack with contents gauge on cylinder regulator |
| HP4021 | 1 | E-size cart system with contents gauge on cylinder regulator |
| HP5102 | 1 | nylon shoulder bag for D cylinder (black) |
| HP5001 | 2 | two-wheel cart for D or E cylinder |
| HP5900 | 1 | plastic wrench for post valve |



Invacare® Lightweight Aluminum Cylinders

Invacare offers empty aluminum oxygen cylinders in standard sizes packaged to meet your requirements. You can order units in packs of six, or in competitively priced pallet quantities.

FEATURES

- All cylinders are available in both CGA870 standard and toggle valves (excluding the CGA540 configured sizes)
- Manufactured by Catalina Cylinders
- Competitively priced
- One-year warranty on post valves

Model Selection Guide

| Cylinder size | Model no. | Units/Case | Oxygen capacity (liters) | Post valve | Valve type |
|---------------|---------------|------------|--------------------------|------------|------------|
| M4 | IRCM4870-6 | 6 | 115 | CGA870 | standard |
| M4 | IRCM4870T-6 | 6 | 115 | CGA870 | toggle |
| M4 | IRCM4870-144 | 144 | 115 | CGA870 | standard |
| M4 | IRCM4870T-144 | 144 | 115 | CGA870 | toggle |
| M6 | IRCM6870-6 | 6 | 164 | CGA870 | standard |
| M6 | IRCM6870T-6 | 6 | 164 | CGA870 | toggle |
| M6 | IRCM6870-144 | 144 | 164 | CGA870 | standard |
| M6 | IRCM6870T-144 | 144 | 164 | CGA870 | toggle |
| ML6 | IRCML6870-6 | 6 | 164 | CGA870 | standard |
| ML6 | IRCML6870T-6 | 6 | 164 | CGA870 | toggle |
| ML6 | IRCML6870-81 | 81 | 164 | CGA870 | standard |
| ML6 | IRCML6870T-81 | 81 | 164 | CGA870 | toggle |
| M7 | IRCM7870-81 | 81 | 198 | CGA870 | standard |
| M7 | IRCM7870T-81 | 81 | 198 | CGA870 | toggle |
| M9 | IRCM9870-6 | 6 | 248 | CGA870 | standard |
| M9 | IRCM9870T-6 | 6 | 248 | CGA870 | toggle |
| M9 | IRCM9870-81 | 81 | 248 | CGA870 | standard |
| M9 | IRCM9870T-81 | 81 | 248 | CGA870 | toggle |
| MD | IRCMD870-6 | 6 | 415 | CGA870 | standard |
| MD | IRCMD870T-6 | 6 | 415 | CGA870 | toggle |
| MD | IRCMD870-81 | 81 | 415 | CGA870 | standard |
| MD | IRCMD870T-81 | 81 | 415 | CGA870 | toggle |
| ME | IRCME870-6 | 6 | 682 | CGA870 | standard |
| ME | IRCME870T-6 | 6 | 682 | CGA870 | toggle |
| ME | IRCME870-81 | 81 | 682 | CGA870 | standard |
| ME | IRCME870T-81 | 81 | 682 | CGA870 | toggle |
| M22 | IRCM22870-72 | 72 | 647 | CGA870 | standard |
| M22 | IRCM22870T-72 | 72 | 647 | CGA870 | toggle |
| M60 | IRCM60540-42 | 42 | 1,726 | CGA540 | handwheel |
| M90 | IRCM90540-42 | 42 | 2,549 | CGA540 | handwheel |
| MM | IRCMM540-25 | 25 | 3,452 | CGA540 | handwheel |



Invacare offers a full line of respiratory products to assist in the treatment of sleep apnea, asthma, chronic obstructive pulmonary disease (COPD) and other respiratory diseases and conditions. The Invacare family of respiratory products is constantly evolving as we strive to provide innovative solutions to meet the ever-changing needs of you and your patients.

We are pleased to share with you our full line Respiratory Catalog. We hope it will become a useful resource for all of your respiratory product needs. From aerosol compressors to sleep apnea products, we've added many significant new products:

- Invacare® Patient Convenient Pack for the Invacare® HomeFill™ II (page 5)
- Invacare® Platinum™ 5 oxygen concentrator (page 6)
- Digit® finger pulse oximeter (page 13)
- Invacare® Poseidon™ passover humidifier (page 16)
- Invacare Freeway® Freedom portable aerosol compressor (page 21)
- VixOne™ nebulizer with Dragon mask or Pedi-Neb™ pacifier (page 23)

We are committed to product innovation as demonstrated in the development of the Invacare HomeFill II oxygen filling system (cover/page 4). This innovative product revolutionizes ambulatory oxygen by allowing patients to fill their own high-pressure cylinders in the convenience of their own home, giving them greater independence and freedom. At the same time, this product was developed to help you benefit from the cost-savings resulting from the virtual elimination of time-consuming and costly service calls associated with cylinder and/or liquid oxygen deliveries.

We share in the concern for the future of the industry and continue to lead the way by lobbying for legislation that supports respiratory providers such as the Advanced Beneficiary Notice.

We are involved in the community and are proud of our sponsorship of the American Association of Respiratory Care (AARC) "Invacare Excellence in Homecare Award," as well as our participation in the American Association of Homecare HME/RT Council and Sleep Disordered Breathing Task Force.