

# Medical Gas Cylinders - Aluminum

## Specifications

Sales Part Number	Service Pressure		Oxygen Capacity		Outside Diameter		Overall Length		Empty Weight		Min. Internal Volume		Min. Water Weight		Thread Size UNF
	psi	bar	cu ft	liters	in	mm	in	mm	lbs	kg	cu in	liters	lbs	kg	
<b>M004</b>	2216	153	4.0	113.3	3.21	81.41	9.0	228.6	1.9	0.8	43	0.7	1.53	0.70	0.750-16 UNF-2B
<b>M006</b>	2216	153	5.8	164.5	3.21	81.41	11.8	299.0	2.2	1.0	62	1.0	2.24	1.01	0.750-16 UNF-2B
<b>M006A</b>	2015	139	6.0	171.1	4.38	111.13	7.9	201.7	2.8	1.2	71	1.2	2.56	1.16	0.750-16 UNF-2B
<b>M007</b>	2015	139	7.4	209.6	4.38	111.13	9.2	233.2	3.3	1.5	87	1.4	3.14	1.42	0.750-16 UNF-2B
<b>M009</b>	2015	139	8.8	248.4	4.38	111.13	10.6	270.0	3.7	1.7	103	1.7	3.71	1.68	0.750-16 UNF-2B
<b>M011</b>	2650	183	11.4	323.7	4.38	111.13	10.9	276.2	4.5	2.0	101	1.7	3.66	1.66	0.750-16 UNF-2B
<b>MD15</b>	2015	139	14.6	414.6	4.38	111.13	16.2	411.2	4.9	2.2	172	2.8	6.20	2.81	0.750-16 UNF-2B
<b>M018</b>	2216	153	18.0	510.9	5.25	133.35	14.1	357.1	7.6	3.4	193	3.2	6.95	3.15	0.750-16 UNF-2B
<b>M022</b>	2216	153	22.9	647.1	5.25	133.35	17.1	435.1	8.9	4.0	244	4.0	8.80	3.99	0.750-16 UNF-2B
<b>ME24</b>	2015	139	24.1	682.2	4.38	111.25	25.2	639.8	7.3	3.3	283	4.6	10.20	4.63	0.750-16 UNF-2B
<b>M060</b>	2216	153	60.9	1723.6	7.25	184.15	23.5	596.9	22.5	10.2	650	10.7	23.44	10.63	1.125-12 UNF-2B, M25 x 2
<b>M00M</b>	2216	153	121.9	3452.5	8.00	203.20	36.3	922.3	40.0	18.1	1302	21.3	46.95	21.30	1.125-12 UNF-2B
<b>M265</b>	2216	153	265.1	7506.8	9.80	248.92	51.9	1319.0	87.0	39.5	2831	46.4	102.08	46.30	0.750-14 NGT, 1.125-12 UNF-2B

3016 Kansas Avenue, Riverside, CA 92507 USA,  
 Tel: (909) 684 5110 or (800) 764 0366 Fax: (909) 781 6598

© 2003 Luxfer Inc.

