



H E A V Y D U T Y S C O O T E R

Seat Color
Charcoal Vinyl

Body Color
Silver

Ride

^{golden}
Avenger
4 wheel

Luxurious, 4 wheel electric scooter features excellent ground clearance, 4 way flashers, oversized seat, and power.

Rise. Roll. Rest. Ride.

GOLDEN
technologies

YOUR LIFE IN MOTION

Avenger 4 wheel Specifications

| | |
|---------------------------------|--|
| Model | GA-541 |
| Length: | 57.5" |
| Width | 24.6" |
| Tires: | |
| Front tire: | 13" |
| Rear Tire: | 13" |
| Seat Depth: | 17.5" |
| Seat Width: | 22" |
| Weight Capacity: | 500 lbs. |
| Maximum Speed: | 7 MPH |
| Ground Clearance: | 4" |
| Operating Range: | up to 18 Miles* |
| Turning Radius: | 67" |
| Batteries (2): | 22 Nf |
| Seat weight: | 58 lbs. |
| Front half weight: | 66 lbs. |
| Rear half weight: | 66 lbs. |
| Freewheel mode: | Yes |
| Electronic speed control: | Yes |
| Brakes: | Electro-mechanical |
| Charger: | Off-board 5 amp |
| Rear wheel drive: | Sealed transaxle, 24 Volt VDC / 4 pole motor |
| Controller: | PG 110 amp / programmable |
| Warranty – Frame: | Limited Lifetime |
| Warranty – Motor/Electronics:** | 2 Years |
| Warranty – Drive Train | 2 Years |



Features

- 500 lb weight capacity
- Speeds up to 7 mph
- Bright safety lights
- 22 Nf batteries
- Up to 18 Mile Range*
- 4-Way Flashers
- Great Warranty
- Oversized tires
- Oversized seat
- Removable basket



* The operational mileage will be less if scooter is equipped with foam filled tires. Battery range will vary based on rider weight and terrain.

** Electronics warranty excludes Batteries.

Literature is current at time of printing. Company reserves the right to make changes.