



F U L L S I Z E S C O O T E R

Seat Color
Black Vinyl

Body Colors
Vermillion Red

Arctic Blue

Ride

^{golden}
COMPANION™ II
3 wheel

Full size, three-wheel scooter offering excellent leg room and comfortable seating.



Rise. Roll. Rest. Ride.



YOUR LIFE IN MOTION

golden COMPANION™ II

GC-322, GC-323, GC-325

F U L L S I Z E S C O O T E R

Companion II 3 wheel Specifications

Model:	GC-322	GC-323	GC-325
Length:	46"	46"	46"
Width:	24"	24"	24"
Tires:	Foam-filled	Foam-filled	Foam-filled
Front Tire:	9"	9"	9"
Rear Tire:	10"	10"	10"
Weight Capacity:	350 lbs.	300 lbs.	350 lbs.
Maximum Speed:	5.5 MPH	5.8 MPH	5.8 MPH
Ground Clearance:	4"	4"	4"
Operating Range:	32 Miles*	32 Miles*	32 Miles*
Turning Radius:	48"	48"	48"
Batteries (2):	Type U-1	Type U-1	Type U-1
Total weight with batteries:	150 lbs.	150 lbs.	150 lbs.
Freewheel mode:	Yes	Yes	Yes
Electronic speed control:	Yes	Yes	Yes
Brakes:	Dynamic/ Electro-mechanical	Dynamic/ Electro-mechanical	Dynamic/ Electro-mechanical
Charger:	Onboard 3.0 amp	Onboard 3.0 amp	Onboard 3.0 amp
Rear wheel drive:	Sealed transaxle, motor 24 Volt	Sealed transaxle, motor 24 Volt	Sealed transaxle, motor 24 Volt
Controller:	70 amp	70 amp	70 amp
Warranty – Frame:	Limited Lifetime	Limited Lifetime	Limited Lifetime
Warranty – Motor/Electronics:**	1 Year	1 Year	1 Year
Warranty – Drive Train:	2 Years	2 Years	2 Years



Patented infinite adjustable tiller



On-board self-diagnostics system

Features

- Easy to Operate
- Terrific Maneuverability
- Great Warranty
- Disassembles for Transportation
- Bright, Metallic Colors

Arctic Blue



Vermillion Red



Disassembles easily! (Four-wheel model shown)

* The operational mileage will be less if scooter is equipped with foam filled tires. Battery range will vary based on rider weight and terrain.

** Electronics warranty excludes Batteries.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

