


F U L L S I Z E S C O O T E R



Seat Color
Black Vinyl

Body Colors
Vermillion Red
Arctic Blue

Ride

golden COMPANION™ II 4 wheel

Full size, 4 wheel scooter offering excellent leg room and comfortable seating.



Rise. Roll. Rest. Ride.



YOUR LIFE IN MOTION

golden COMPANION™ II

GC-422

F U L L S I Z E S C O O T E R

Companion II 4 wheel Specifications

Model	GC-422
Length:	51"
Width	23"
Tires:	Foam-filled
Front tire:	Mag wheel 9"
Rear Tire:	10"
Weight Capacity:	350 lbs.
Maximum Speed:	5.9 MPH
Ground Clearance:	4"
Operating Range:	up to 19 Miles*
Turning Radius:	54"
Batteries (2):	Type U-1
Total weight	with batteries 179 lbs.
Freewheel mode:	Yes
Electronic speed control:	Yes
Brakes:	Dynamic/Electro-mechanical
Charger:	Onboard 3.0 amp
Rear wheel drive:	Sealed transaxle, motor 24 Volt
Controller:	70 amp
Warranty – Frame:	Limited Lifetime
Warranty – Motor/Electronics:**	1 Year
Warranty – Drive Train	2 Years



Patented infinite adjustable tiller



On-board self-diagnostics system

Features

- Disassembles for transportation
- Bright, Metallic Colors



Easy to operate



Independent front suspension



Disassembles easily!

* The operational mileage will be less if scooter is equipped with foam filled tires. Battery range will vary based on rider weight and terrain.

** Electronics warranty excludes Batteries.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

